



Once an alarm has been triggered the alarm capsule with the balloon **cannot be reused** again. Place a new alarm capsule into the monitoring bracelet. **Attention:** capsule needs also to be replaced when reaching the expiry date marked on the capsule. **Before inserting a new capsule, warm the holder (without capsule) in the oven at 70 degrees for approx. 10 min.** ① Make sure that the balloon is almost flat, before placing it into the bracelet. There might be some air left in the balloon. Lightly squeeze the balloon for air release. ② The docking part with the pins on the capsule has to be placed at the LED part of the bracelet. On the **capsule** the arrow has to be placed direction LED window. The field with the expiry date, for example 03/2020 has to point in direction of the short bracelet part with the buckle and pin. **Attention:** Twist the capsule slightly and use little pressure in order to bring the capsule into the bracelet until the green LED appears. ③ The green LED will blink each 5 seconds and now BlueFox ST1 is ready to be used.